

RK-06F



Brand:
Product Code: RK-06F

Short Description

Portable rack
Fast assembly
Fine-focus adjustment

Description

The Dino-Lite RK-06F stand is a sturdy and stable mid-range stand solution, with fine focus adjustment that can be used with all Dino-Lite digital microscopes. It is constructed of stainless steel and lightweight aluminum and offers precise fine-focus adjustment as well as a quick vertical release function.

Its adjustable horizontal arm has a 15cm range with 360-degree rotation and with the RK-10-EX arm (excluded) its effective range can be extended to observe larger objects. Various optional accessories are available for the RK-06F, such as the RK-10-FX flexible arm or the RK-10-PX, XY arm. With these features, the RK-06F is the ideal mid-range stand solution for precise orientation and maximum flexibility.

[AVAILABLE EXTENSIONS FOR THE RK-06A STAND](#)

RK-10-EX



Additional horizontal arm (included with RK-10A)

RK-10-FX



Flexible arm extension

RK-10-VX



Vertical arm extension

RK-10-PX2



X and Y positioning arm

Specification

Compatibility	
Compatible with:	All Dino-Lite models
Housing	
Dimensions	22 cm (L) x 15 cm (W) x 27 cm (H)
Weight	1060 g
Features	
ESD safe	No

Information	
Package contents	Stand, HD-P1 holder, RK-10-EX horizontal arm
Warranty information	2 years European warranty
Accessory range:	Professional stands
Price range	€250,00 - €400,00

Product Gallery

