

RK-05F

Brand: Dino-Lite
Product Code: RK-05F



Short Description

Portable rack
Tilt and stay feature (30 degrees)
Fast assembly
Fine-focus adjustment

Description

RK-05F is a light-weight, ultra-portable stand with a one-step assemble/disassemble design.

If you want to observe your object or specimen from a different perspective, the Dino-Lite RK-05F stand provides a multi-angle view with its 30 degrees inclinable pole feature. Additionally it has an integrated fine focus adjustment knob to make focusing easier and more precise.

RK-05F is a light-weight, ultra-portable stand with a one-step assemble/disassemble design.

All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-05 without fine focus adjustment](#)

Specification

Compatibility	
Compatible with:	All Dino-Lite models, except Long Working Distance models
Housing	
Dimensions	14cm (W) x 11cm (L) x 16.5cm (H)
Weight	636 g
Features	
ESD safe	No
Information	
Package contents	Stand, HD-P1 holder
Warranty information	2 years European warranty
Accessory range:	Professional stands
Price range	€200,00 - €300,00

Product Gallery

