
Dino-Lite Europe - Phone: - Email:

RK-10F / RK-10FS





Short Description

Quick Release Button Fine-Focus Knob Anti-Backlash Lock Arm Screw HD-M2A

Description

The RK-10F is an upgraded version of the RK-10, featuring an improved quick release button with an anti-backlash mechanism that enhances its ability to make precise focus adjustments, ensuring smoother and more responsive operation.

The RK-10FS incorporates a heavier antistatic baseplate that suppresses external vibrations, making it an optimal choice for stable operation under high-magnification inspections.

AVAILABLE EXTENSIONS FOR THE RK-10A STAND

RK-10-EX



Additional horizontal arm (included with RK-10A)

RK-10-FX



Flexible arm extension

RK-10-VX



Vertical arm extension

RK-10-PX2



X and Y positioning arm

Specification

Housing	
Dimensions	226 mm (H) x 150 mm (D) x 220 mm (L)
Dimensions	226 mm (H) x 150 mm (D) x 220 mm (L)
Weight	Approx. 1050 g
Weight	Approx. 1900 g
Features	
ESD safe	Yes

Information	
Package contents	Stand, HD-P1 holder, RK-10-EX horizontal arm
Warranty information	2 years European warranty
Accessory range:	Professional stands
Price range	€300,00 - €450,00

Product Gallery

