

SW-F1

Brand: Dino-Lite
Product Code: SW-F1



Short Description

Footpedal

Description

The SW-F1 allows you to take a picture with a simple step from your foot using the Dino-Lite or Dino-Eye digital microscopy solutions. The foot pedal is only compatible with DinoCapture 2.0 for the PC and DinoXcope for Macintosh computers. No additional software is required, simply plug in the device into your USB port and it is ready to use.

Specification

Compatibility	
Compatible with:	All USB models, except Wireless models
Housing	
Dimensions	10cm (W) x 7cm (L) x 3cm (H)
Weight	300 g
Features	
ESD safe	No
Information	
Package contents	1 foot pedal
Warranty information	2 years European warranty
Accessory range:	Light and control
Price range	€50,00 - €150,00

Product Gallery

