

# RK-06A



**Brand:**  
**Product Code:** RK-06A

## Short Description

Portable rack  
Fast assembly  
Fine-focus adjustment

## Description

The Dino-Lite RK-06A stand is a sturdy and stable mid-range stand solution that can be used with all Dino-Lite digital microscopes. It is constructed of stainless steel and lightweight aluminum and offers precise fine-focus adjustment as well as a quick vertical release function.

Its adjustable horizontal arm has a 15cm range with 360-degree rotation and with the RK-10-EX arm (included) its effective range can be extended to observe larger objects. Various optional accessories are available for the RK-06A, such as the RK-10-FX flexible arm or the RK-10-PX, XY arm. With these features, the RK-06A is the ideal mid-range stand solution for precise orientation and maximum flexibility.

[AVAILABLE EXTENSIONS FOR THE RK-06A STAND](#)

---

**RK-10-EX**



Additional horizontal arm (included with RK-10A)

**RK-10-FX**



Flexible arm extension

**RK-10-VX**



Vertical arm extension

**RK-10-PX2**



X and Y positioning arm

### Specification

|                  |                                   |
|------------------|-----------------------------------|
| Compatibility    |                                   |
| Compatible with: | All Dino-Lite models              |
| Housing          |                                   |
| Dimensions       | 22 cm (L) x 15 cm (W) x 27 cm (H) |
| Weight           | 1060 g                            |
| Features         |                                   |
| ESD safe         | No                                |

## Information

|                      |  |
|----------------------|--|
| Package contents     | Stand, HD-P1 holder, RK-10-EX horizontal arm |
| Warranty information | 2 years European warranty                    |
| Accessory range:     | Professional stands                          |
| Price range          | €200,00 - €300,00                            |

## Product Gallery

