

# KM-01

**Brand:** Dino-Lite  
**Product Code:** KM-01



## Short Description

Remote focus control

## Description

The KM-01 can be used to control the Dino-Lite focus wheel without touching the device. This is especially useful when the Dino-Lite needs to be isolated during use.

The touch controls on the durable KM-01 control box provide instant response, but the KM-01 can also be controlled with the DinoCapture 2.0 software. The KM-01 can be mounted on the RK-04/RK-05/RK-06A/RK-10A Dino-Lite stands, or on 3rd party stands with a 10mm diameter pole.

*The KM-01 is compatible with certain models of Dino-Lite. Please consult your local reseller or Dino-Lite Europe for further information.*

## Specification

Compatibility	
Compatible with:	Please inquire with your local Dino-Lite sales representative for compatibility information
Housing	
Dimensions	77.14mm (L) 65.55mm (W) 98mm (H)
Weight	340 g (Motor) 110 g (Control Box)

Features	
ESD safe	No
Information	
Package contents	Motor, control box, c-clip, set screws, AC adapter, USB cable, male/female 2 pin circular DIN connector, software CD
Warranty information	2 years European warranty
Accessory range:	Light and control
Price range	€400,00 - €600,00

## Product Gallery

