

MS53BA2

Brand: Dino-Lite
Product Code: MS53BA2

Short Description

Flexible arm stand with table base

Description

Heavy duty table base for the MS53B. The MS53BA2 is a combination of the rigid MS53B articulating flex-arm and a desktop base rack that holds the MS53B on top of a table.

The desktop base rack is designed for situations where the standard clamp cannot be used to secure the MS53B on a stable surface. The flex-arm is moved and adjusted with tightening knobs and clamps for positioning the arm to the area you wish to examine. Two of the flex-arm axis can be rotated 360 degrees and extended for further reach. The MS53BA2 flex-arm is manufactured from durable, lightweight recycled aluminum and has the ability to be folded flat, which makes the MS523BA2 stand portable and ready for travel.

Specification

| | |
|----------------------|-----------------------------------|
| Compatibility | |
| Compatible with: | All Dino-Lite models |
| Housing | |
| Dimensions | 35.5cm (W) x 16.5cm (L) x 5cm (H) |
| Weight | 4 kg |
| Features | |
| ESD safe | No |
| Information | |
| Package contents | Stand, table base, HL2 holder |
| Warranty information | 2 years European warranty |
| Accessory range: | Professional stands |
| Price range | €400,00 - €600,00 |

Product Gallery

