

F25MM

Brand: Dino-Lite
Product Code: F25MM



Short Description

C-mount adapter for rigid endoscope

Description

C-mount adapter for rigid endoscopes with a focal length of 25mm. This C-mount adapter connects to rigid, straight through endoscopes with max. 32mm ocular diameter.

Specification

Compatibility	
Compatible with:	AM4023CT, AM4025X, AM7025X
Housing	
Dimensions	48mm (H) x 52mm (W)
Weight	120 g
Features	
ESD safe	No
Information	
Package contents	Adapter
Warranty information	2 years European warranty
Accessory range:	Various
Price range	€250,00 - €350,00

Product Gallery

