Dino-Lite Europe - Phone: - Email:

CS-41

Product Code: CS-41



Short Description

Glass calibration slide

Description

The CS-41 is a calibration slide that can be used to calibrate your Dino-Lite Edge digital microscope or to be used with the Dino-Eye Edge eyepiece cameras. The total length is 10mm (10.000um) and the smallest gradation on the slide is 0.1mm (100um). The CS-41 calibration slide also contains QR codes in 3 different sizes that can be used for the Auto Calibration procedure within DinoCapture 2.0 software version 1.5.25 or newer. Auto calibration provides an easy way to calibrate a Dino-Lite Edge serie microscope. The DinoCapture software identifies the special QR codes and performs the calibration automatically.

Row A: magnifications of 120x and above

Row B: magnifications of 40x - 220x Row C: magnifications of 20x - 50x

Click here to download the Auto Calibration procedure

Specification

Compatibility	
Compatible with:	All Dino-Lite models

Housing	
Dimensions	75mm (L) x 25mm (W) x 1,5mm (H)
Weight	10 g
Features	
ESD safe	No
Information	
Package contents	1 piece
Warranty information	2 years European warranty
Accessory range:	Various
Price range	€0 - €100,00