

# SL-ZW1

**Brand:** Dino-Lite  
**Product Code:** SL-ZW1



## Short Description

Polarized Sidelight

## Description

The SL-ZW1 is a polarized USB powered gooseneck light with the ability for the polarization to be adjusted by a simple rotation on the accessory itself. To take advantage of the polarizing light, the accessory is ideally suited to work with Dino-Lite models that have the polarization filter.

## Specification

Compatibility	
Compatible with:	All Dino-Lite models
Housing	
Dimensions	40cm (L) x 1.6cm (D)
Weight	56 g
Features	
ESD safe	No
Information	
Package contents	1 side light
Warranty information	2 years European warranty
Accessory range:	Light and control
Price range	€0 - €100,00

## Product Gallery

