

# FC-L-MA1



**Brand:** Dino-Lite  
**Product Code:** FC-L-MA1

## Short Description

Mirror angle adapter

## Description

The mirror angle adapter fits onto the models with long working distance:

AD4113TL/AD4116TL/AD4112PTL/AD4013MTL/AD4013MZTL/AD7013MTL etc.

Mirror angle: 45 °

## Specification

Compatibility	
Compatible with:	All AD models with long working distance
Housing	
Dimensions	68mm (L) x 7mm (W/o tip)
Weight	9 g
Features	
ESD safe	Yes
Information	
Package contents	1 cap
Warranty information	2 years European warranty
Accessory range:	Various
Price range	€50,00 - €150,00

Product Gallery

