

PH-DESKSET – 10044

Brand: Photonic
Product Code: 10044

Short Description

High power LED spot, 2 pieces
 Flexible arms, 2 individual pieces
 Individually mounted on a base plate
 Control unit: brightness and on/off

Description

The Photonic Desk-set is composed of the LED Hi-Power Spot (2 pieces), 2 spot holders (flexible arms) and 2 base plates. The Hi-Power Spot must be attached to the included Spot holder. This Spot holder is then mounted onto the round Base plate. With this set you can create two individual Hi-power spots with flexible arms. As it is mounted on the baseplate, it can be placed anywhere around the sample creating total flexibility. The brightness can be controlled by the compact control unit. Furthermore you can switch on/off the 2-arms individually.

** All other Photonic items are also available through Dino-Lite Europe.*

Specification

Compatibility	
Compatible with:	All Dino-Lite models, except high magnification models
Housing	
Dimensions	55cm long, each arm
Weight	125 g
Features	
ESD safe	Yes
Information	
Package contents	Hi-Power spot arm (2x), spot holders (2x), base plate (2x), controller unit, power supply with EU & UK plug ESD safe add-ons not included
Warranty information	2 years European warranty
Accessory range:	Light and control

Price range	€450,00 - €550,00
-------------	-------------------

Product Gallery

