

WF-10

Brand: Dino-Lite
Product Code: WF-10

Short Description

Dino-Lite Wi-Fi streamer
 Live on mobile devices
 Camera settings
 Multi-users
 Stationary use

Description

The WF-10 WiFi Streamer allows Dino-Lite USB models to be used wirelessly in connection with an iOS (iPhone/iPad) or Android app on any tablet, smartphone, or computer. The WiFi Streamer is ideal for field work or presentations. Stream from one Dino-Lite to multiple devices. The WF-10 has a replaceable battery and can be used anywhere to create a WiFi signal for streaming the Dino-Lite live images.

*(*The WF-10 is compatible with certain models of Dino-Lite. Please consult your local reseller or Dino-Lite Europe for further information.)*

[BC-120B: Spare battery and charger for WF-10.](#)



Specification

| Compatibility | |
|------------------|---|
| Operating system | Windows 7, 8, 10 & 11, iOS 8.0 or later, Android 4.0 or later |
| Compatible with: | AM/AD4xx3, AM/AD7xx3, MEDL4xx, MEDL7Xxx all (R4)-series or newer, AM4xx5, AM7xx5 series * |
| App | DinoConnect |

| Housing | |
|----------------------|---|
| Dimensions | 12.1cm (W) x 5.5cm (L) x 1.6cm (H) |
| Weight | 112 g |
| Features | |
| ESD safe | No |
| Information | |
| Package contents | WF-10 WiFi streamer, user manual, power adapter, 1880mAh Li-Ion battery |
| Warranty information | 2 years European warranty |
| Accessory range: | Wireless |
| Price range | €250,00 - €400,00 |
| Note | *USB 3.0 and EdgePLUS series are not compatible. |

Product Gallery

