

RK-04



Brand: Dino-Lite
Product Code: RK-04

Short Description

Portable rack
Fast assembly

Description

The Dino-Lite RK-04 stand has a robust, lightweight and compact rack design. It is designed with portability and convenience in mind. All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-04F with fine focus adjustment](#)

Specification

Compatibility	
Compatible with:	All Dino-Lite models, except Long Working Distance models
Housing	
Dimensions	14cm (W) x 11cm (L) x 16.5cm (H)
Weight	605 g
Features	

ESD safe	No
Information	
Package contents	Stand, HD-P1 holder
Warranty information	2 years European warranty
Accessory range:	Professional stands
Price range	€100,00 - €200,00

Product Gallery

