

# EarScope Pneumatic (MEDL4EP)

**Brand:** Dino-Lite Medical  
**Product Code:** MEDL4EP



## Short Description

1.3 Megapixel resolution  
Pro model for pneumatic otoscopy  
55-90x Magnification  
Medical Device Class 1

## Description

The EarScope Pneumatic (MEDL4EP) is a compact and safe solution for watching the eardrum and the external ear. The MEDL4EP has a 1,3 Megapixel resolution and a magnification of 55 ~ 90x. It comes with a rubber bulb that is squeezed to give a puff of air into the ear canal, which allows the doctor to see how the eardrum moves. This model is delivered with disposable specula of 2,5 and 4 mm (10 each).



USB 2.0



6 White LEDs  
switchable



Adjustable  
~55x-90x



1.3 megapixel  
1280 x 1024



Magnification  
lock

## Specification

Lighting

Light/ LED type	White
Number of LEDs	6
LED on/off switchable:	Yes
Infrared filter	IR cut-filter >650 nm
Diffuser available	No
Polarizer	No
<b>Optics</b>	
Magnification	55-90x
Lens type	Glass with anti-reflection coating
<b>Sensor</b>	
Sensor type	CMOS
Resolution	1.3 Megapixel (1280x1024)
Maximum frame rate	30 fps
<b>Compatibility</b>	
Interface	USB 2.0
Operating system	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Included software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Supported image formats (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Supported video formats (Windows)	WMV, FLV, SWF
Supported image formats (MacOS)	JPEG, PNG
Supported video formats (MacOS)	MOV
Imaging standards	DirectShow, UVC
<b>Housing</b>	
Housing material	Composite/plastic housing
Magnification lock	Yes
Dimensions	10.5cm (L) x 3.2cm (D)
Weight	95 g
Cable length	1.8m
<b>Features</b>	
Special feature	No
Measurement	No
Calibration	No
Microtouch sensor	Yes
<b>Information</b>	
Package contents	Microscope, carry pouch, 10x 2,5 mm specula, 10x 4 mm specula, software CD, user manual, rubber bulb

Warranty information	2 years European warranty
Regulatory approval	Medical Device Class 1 – Medical Devices Regulation (EU) 2017/745
Price range	€600,00 - €800,00

### Product Gallery

