

# MSAA502



**Brand:** Dino-Lite  
**Product Code:** MSAA502

## Short Description

Multipurpose adapter/ light shield

## Description

Useful for eliminating ambient light or providing an enhanced field of view at optimal distances for various applications. The MSAA502 is not compatible with polarizer models or models in the AD range.

## Specification

Compatibility	
Compatible with:	Dino-Lite AM models, excluding models with polarizer
Housing	
Dimensions	124mm (H) x 84mm (D)
Weight	100 g
Features	
ESD safe	No
Information	
Package contents	1 multipurpose adapter
Warranty information	2 years European warranty
Accessory range:	Various
Price range	€0 - €100,00

Product Gallery



