

# MSAA502



**Brand:** Dino-Lite  
**Product Code:** MSAA502

## Short Description

Multipurpose adapter/ light shield

## Description

Useful for eliminating ambient light or providing an enhanced field of view at optimal distances for various applications. The MSAA502 is not compatible with polarizer models or models in the AD range.

## Specification

<b>Compatibility</b>	
Compatible with:	Dino-Lite AM models, excluding models with polarizer
<b>Housing</b>	
Dimensions	124mm (H) x 84mm (D)
Weight	100 g
<b>Features</b>	
ESD safe	No
<b>Information</b>	
Package contents	1 multipurpose adapter
Warranty information	2 years European warranty
Accessory range:	Various
Price range	€0 - €100,00

# Product Gallery



