

N3C-D Diffuser Cap

Brand: Dino-Lite
Product Code: N3C-D



Short Description

Usage: Working with shiny objects (metal/plastics/PCB/jewelry etc.).

Description

The interchangeable front caps are designed to be used with Dino-Lite Edge and Edge^{PLUS} series models. The caps provide additional versatility to Dino-Lite digital microscope applications. This is a translucent cap that partially softens the LED light. When working with shiny objects (metal/plastics/PCB/jewelry, etc)

Compatibility

All Dino-Lite Edge/Edge^{PLUS} models

N3C caps overview

	MODEL	VERSIO N	N3C-O	N3C-E	N3C-C	N3C-L	N3C- D/D2	N3C-S
			Open cap	Extended open cap	Closed cap	Long cap	Diffuser cap	Sidelight cap
AM	4115, 4117, 4515, 4517, 4815, 491	Standard models: ZT, ZTL, MZT, MZTL,	?	?	?	?	?	?

	7, 7115, 7515, 7915, 73115, 73515, 73915	ZTW						
AM	4115, 4515, 4117, 4517, 7115	Ultraviolet /Infrared models: FVW, FUW, FIT, FJT, FKT , FUT, FVT,	? (black version)	x	x	x	x	x
AM	4115, 4117, 4515, 4517, 4815	Fluorescen t models: CFVW, DFRW, GFBW, G2FBW, GRFBY, YFGW, RFYW	? (black version)	x	x	x	x	x
AM	4515, 7515, 73515	High magn ification models: ZT4, T5, T8, MT4A, MT8A	?	x	?	x	x	x
AM	4115, 7115, 73115	Extra Long Working distance models: ZT, ZTF	?	x	x	x	x	x
AM	5216, 5218	VGA, HDMI models: ZT, ZTL	?	?	?	?	?	?