

N3C-O Open cap

Brand: Dino-Lite
Product Code: N3C-O



Short Description

Usage: Generic use

Description

The interchangeable front caps are designed to be used with Dino-Lite Edge and Edge ^{PLUS} series models. The caps provide additional versatility to Dino-Lite digital microscope applications. This 9.5mm cap eases the task of finding and holding a fixed focus when using the Dino-Lite in handheld mode.

Compatibility

All Dino-Lite Edge/ Edge ^{PLUS} models

N3C caps overview

	MODEL	VERSIO N	N3C-O	N3C-E	N3C-C	N3C-L	N3C- D/D2	N3C-S
			Open cap	Extended open cap	Closed cap	Long cap	Diffuser cap	Sidelight cap
AM	4115, 4117, 4515, 4517, 4815, 491 7, 7115,	Standard models: ZT, ZTL, MZT, MZTL, ZTW	?	?	?	?	?	?

	7515, 7915, 73115, 73515, 73915							
AM	4115, 4515, 4117, 4517, 7115	Ultraviolet /Infrared models: FVW, FUW, FIT, FJT, FKT , FUT, FVT,	? (black version)	x	x	x	x	x
AM	4115, 4117, 4515, 4517, 4815	Fluorescen t models: CFVW, DFRW, GFBW, G2FBW, GRFBY, YFGW, RFYW	? (black version)	x	x	x	x	x
AM	4515, 7515, 73515	High magn ification models: ZT4, T5, T8, MT4A, MT8A	?	x	?	x	x	x
AM	4115, 7115, 73115	Extra Long Working distance models: ZT, ZTF	?	x	x	x	x	x
AM	5216, 5218	VGA, HDMI models: ZT, ZTL	?	?	?	?	?	?