

# MS12C

**Brand:** Dino-Lite

**Product Code:** MS12C



## Short Description

Transparent holder

## Description

Variable view range caddy with two variable view range adapters, transparent. The MS12C is not compatible with polarizer models or models in the AD Pro 2 range. The MS12C consists of a base unit with two insert adapters providing three "hands free" fixed working distance and magnification options.

Ideal for collectors, hobbyists, graphic artists as well as at home or in the classroom. The Dino-Lite digital microscopes simply slide into the base or one of the adapters. The clear acrylic base design allows easy accessibility of tools and light to the observed object.

## Specification

Compatibility	
Compatible with:	All Dino-Lite models, except models AD-model, Long Working Distance models and models with polarizer
Housing	
Dimensions	8.2cm (W) x 8.2cm (L) x 7.6cm (H)
Weight	54 g
Features	

ESD safe	No
Information	
Package contents	Stand + 2 adapters
Warranty information	2 years European warranty
Accessory range:	Basic stands
Price range	€0 - €100,00

## Product Gallery

