

PH-HPS SET

Brand: Photonic

Product Code: 10068 + 619-20-049

Short Description

High power LED spots

2-arms, flexible

RK-06A/RK-06AE/RK-06F stand adapter

Control unit: brightness and on/off

Description

The Photonic HPS-set is composed of the LED Hi-Power Spot 2-arm and the RK-06A/RK-06AE/RK-06F adapters.

This spot is used for direct attachment to microscope stands, like the Dino-Lite RK-06A/RK-06AE/RK-06F stands. As it is mounted on the stand, it is perfect for use in small spaces. The brightness can be controlled by the compact control unit. Furthermore you can switch on/off the 2-arms individually. The flexible arms can be easily adjusted to provide the best angle of illumination to your sample.

** All other Photonic items are also available through Dino-Lite Europe.*

Specification

Compatibility	
Compatible with:	All Dino-Lite models, except high magnification models
Housing	
Dimensions	55cm long, each arm
Weight	650 g
Features	
ESD safe	Yes
Information	
Package contents	Hi-Power Spot 2-arms (10068) , controller unit, power supply with EU & UK plug, RK-06A/RK-06AE/RK-06F adapter (619-20-049) ESD Safe add-ons not included

Warranty information	2 years European warranty
Accessory range:	Light and control
Price range	€450,00 - €550,00

Product Gallery

