

RK-04F



Brand: Dino-Lite
Product Code: RK-04F

Short Description

Portable rack
Fast assembly
Fine-focus adjustment

Description

The Dino-Lite RK-04F stand has a robust, lightweight and compact rack design. Additionally it has an integrated fine focus adjustment knob to make focusing easier and more precise. The RK-04F is the perfect accessory for users requiring flexibility.

All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-04 without fine focus adjustment](#)

Specification

Compatibility	
Compatible with:	All Dino-Lite models, except Long Working Distance models
Housing	

Dimensions	14cm (W) x 11cm (L) x 16.5cm (H)
Weight	630 g
Features	
ESD safe	No
Information	
Package contents	Stand, HD-P1 holder
Warranty information	2 years European warranty
Accessory range:	Professional stands
Price range	€150,00 - €250,00

Product Gallery

