

PH-DESKSET – 10044

Brand: Photonic
Product Code: 10044

Short Description

High power LED spot, 2 pieces
Flexible arms, 2 individual pieces
Individually mounted on a base plate
Control unit: brightness and on/off

Description

The Photonic Desk-set is composed of the LED Hi-Power Spot (2 pieces), 2 spot holders (flexible arms) and 2 base plates. The Hi-Power Spot must be attached to the included Spot holder. This Spot holder is then mounted onto the round Base plate. With this set you can create two individual Hi-power spots with flexible arms. As it is mounted on the baseplate, it can be placed anywhere around the sample creating total flexibility. The brightness can be controlled by the compact control unit. Furthermore you can switch on/off the 2-arms individually.

** All other Photonic items are also available through Dino-Lite Europe.*

Specification

| | |
|----------------------|---|
| Compatibility | |
| Compatible with: | All Dino-Lite models, except high magnification models |
| Housing | |
| Dimensions | 55cm long, each arm |
| Weight | 125 g |
| Features | |
| ESD safe | Yes |
| Information | |
| Package contents | Hi-Power spot arm (2x), spot holders (2x), base plate (2x), controller unit, power supply with EU & UK plug ESD safe add-ons not included |
| Warranty information | 2 years European warranty |
| Accessory range: | Light and control |

Price range

€450,00 - €550,00

Product Gallery

