

RK-04



Značka: Dino-Lite
Kód produktu: RK-04

Short Description

Portable rack
Fast assembly

Popis

The Dino-Lite RK-04 stand has a robust, lightweight and compact rack design. It is designed with portability and convenience in mind. All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-04F with fine focus adjustment](#)

Specifikace

| | |
|------------------|---|
| Compatibility | |
| Compatible with: | All Dino-Lite models, except Long Working Distance models |
| Housing | |
| Dimensions | 14cm (W) x 11cm (L) x 16.5cm (H) |
| Weight | 605 g |
| Features | |

| | |
|----------------------|---------------------------|
| ESD safe | No |
| Information | |
| Package contents | Stand, HD-P1 holder |
| Warranty information | 2 years European warranty |
| Accessory range: | Professional stands |
| Price range | €100,00 - €200,00 |

Product Gallery

