

CA1070

Brand: Dino-Lite
Product Code: CA1070



Short Description

Sturdy watertight carrying case

Description

Sturdy, watertight case for Dino-Lite microscopes. This heavy duty carrying case with foam inlay protects your Dino-Lite against dust, scratches and moisture. Watertight according to IP67 standard.

Specification

| | |
|----------------------|---|
| Compatibility | |
| Compatible with: | All Dino-Lite models |
| Housing | |
| Dimensions | Exterior: 210 x 167 x 90mm Interior: 186 x 123 x 75mm |
| Weight | 570 g |
| Features | |
| ESD safe | No |
| Information | |
| Package contents | Case, foam inlay |
| Warranty information | 2 years European warranty |
| Accessory range: | Various |
| Price range | €0 - €50,00 |

Product Gallery

