

RK-06A



Marke:
Artikelnummer: RK-06A

Short Description

Portable rack
Fast assembly
Fine-focus adjustment

Beschreibung

The Dino-Lite RK-06A stand is a sturdy and stable mid-range stand solution that can be used with all Dino-Lite digital microscopes. It is constructed of stainless steel and lightweight aluminum and offers precise fine-focus adjustment as well as a quick vertical release function.

Its adjustable horizontal arm has a 15cm range with 360-degree rotation and with the RK-10-EX arm (included) its effective range can be extended to observe larger objects. Various optional accessories are available for the RK-06A, such as the RK-10-FX flexible arm or the RK-10-PX, XY arm. With these features, the RK-06A is the ideal mid-range stand solution for precise orientation and maximum flexibility.

[AVAILABLE EXTENSIONS FOR THE RK-06A STAND](#)

RK-10-EX



Additional horizontal arm (included with RK-10A)

RK-10-FX



Flexible arm extension

RK-10-VX



Vertical arm extension

RK-10-PX2



X and Y positioning arm

Spezifikation

| | |
|------------------|-----------------------------------|
| Compatibility | |
| Compatible with: | All Dino-Lite models |
| Housing | |
| Dimensions | 22 cm (L) x 15 cm (W) x 27 cm (H) |
| Weight | 1060 g |
| Features | |
| ESD safe | No |

Information

| | |
|----------------------|--|
| Package contents | Stand, HD-P1 holder, RK-10-EX horizontal arm |
| Warranty information | 2 years European warranty |
| Accessory range: | Professional stands |
| Price range | €200,00 - €300,00 |

Product Gallery

