

# IriScope Wireless (MEDLW4R)



**Marke:** Dino-Lite Medical  
**Artikelnummer:** MEDLW4R

## Short Description

1.3 Megapixel resolution  
Make detailed pictures of the Iris  
10-20x Magnification  
Wireless  
Medical Device Class 1

## Beschreibung

The IriScope Wireless (MEDLW4R) offers a new, fast and easy way to make pictures of eyes. The IriScope Wireless has 2x2 LED's (white/yellow) that illuminate the side of the eye which reduces the reflections in the image as well as the discomfort by bright light. The yellow LED's provide better images of dark irises. The images have a maximum resolution of 1.3 Megapixel.



USB 2.0



Switchable  
white/yellow  
LEDs



Adjustable  
~10-20x



Magnification  
lock



1.3 megapixel  
1280 x 1024



Wireless

## Spezifikation

Lighting

Light/ LED type	White/ Yellow
Number of LEDs	2+2
LED on/off switchable:	Yes
Infrared filter	IR cut-filter >650 nm
Diffuser available	No
Polarizer	No
<b>Optics</b>	
Magnification	10-20x
Lens type	Glass with anti-reflection coating
<b>Sensor</b>	
Sensor type	CMOS
Resolution	1.3 Megapixel (1280x1024)
Maximum frame rate	30 fps
<b>Compatibility</b>	
Interface	USB 2.0
Operating system	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Supported image formats (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Supported video formats (Windows)	WMV, FLV, SWF
Supported image formats (MacOS)	JPEG, PNG
Supported video formats (MacOS)	MOV
Imaging standards	DirectShow, UVC
Wifi	Yes
<b>Housing</b>	
Housing material	Composite/plastic housing
Magnification lock	Yes
Dimensions	10.5cm (L) x 3.2cm (D)
Weight	95 g
Cable length	1.8m
<b>Features</b>	
Special feature	No
Measurement	No
Calibration	No
Microtouch sensor	Yes
<b>Information</b>	
Package contents	Microscope, carry pouch, software CD, user manual

Warranty information	2 years European warranty
Regulatory approval	Medical Device Class 1 – Medical Devices Regulation (EU) 2017/745
Price range	€900,00 - €1100,00

## Product Gallery

