

# MS08B

**Marke:** Dino-Lite  
**Artikelnummer:** MS08B



## Short Description

Compact basic stand

## Beschreibung

The MS08B is an ultra-portable mini stand that holds the Dino-Lite in place whenever you need to. The MS08B is a great little accessory to enhance the Dino-Lite's microscopic portability solution.

## Spezifikation

<b>Compatibility</b>	
Compatible with:	All Dino-Lite models, except Extra Long Working Distance models
<b>Housing</b>	
Dimensions	7.5cm (W) x 10cm (L) x 14,5cm (H)
Weight	500 g
<b>Features</b>	
ESD safe	No
<b>Information</b>	
Package contents	Stand, HD-P1 holder
Warranty information	2 years European warranty
Accessory range:	Basic stands
Price range	€50,00 - €150,00

**Product Gallery**

