

MS08B

Marke: Dino-Lite
Artikelnummer: MS08B



Short Description

Compact basic stand

Beschreibung

The MS08B is an ultra-portable mini stand that holds the Dino-Lite in place whenever you need to. The MS08B is a great little accessory to enhance the Dino-Lite's microscopic portability solution.

Spezifikation

Compatibility	
Compatible with:	All Dino-Lite models, except Extra Long Working Distance models
Housing	
Dimensions	7.5cm (W) x 10cm (L) x 14,5cm (H)
Weight	500 g
Features	
ESD safe	No
Information	
Package contents	Stand, HD-P1 holder
Warranty information	2 years European warranty
Accessory range:	Basic stands
Price range	€50,00 - €150,00

Product Gallery

