

# MS53BA2

**Marke:** Dino-Lite

**Artikelnummer:** MS53BA2

## Short Description

Flexible arm stand with table base

## Beschreibung

Heavy duty table base for the MS53B. The MS53BA2 is a combination of the rigid MS53B articulating flex-arm and a desktop base rack that holds the MS53B on top of a table.

The desktop base rack is designed for situations where the standard clamp cannot be used to secure the MS53B on a stable surface. The flex-arm is moved and adjusted with tightening knobs and clamps for positioning the arm to the area you wish to examine. Two of the flex-arm axis can be rotated 360 degrees and extended for further reach. The MS53BA2 flex-arm is manufactured from durable, lightweight recycled aluminum and has the ability to be folded flat, which makes the MS523BA2 stand portable and ready for travel.

## Spezifikation

<b>Compatibility</b>	
Compatible with:	All Dino-Lite models
<b>Housing</b>	
Dimensions	35.5cm (W) x 16.5cm (L) x 5cm (H)
Weight	4 kg
<b>Features</b>	
ESD safe	No
<b>Information</b>	
Package contents	Stand, table base, HL2 holder
Warranty information	2 years European warranty
Accessory range:	Professional stands
Price range	€400,00 - €600,00

## Product Gallery

