

# MS53BA2

**Marca:** Dino-Lite

**Código del producto:** MS53BA2

## Short Description

Flexible arm stand with table base

## Descripción

Heavy duty table base for the MS53B. The MS53BA2 is a combination of the rigid MS53B articulating flex-arm and a desktop base rack that holds the MS53B on top of a table.

The desktop base rack is designed for situations where the standard clamp cannot be used to secure the MS53B on a stable surface. The flex-arm is moved and adjusted with tightening knobs and clamps for positioning the arm to the area you wish to examine. Two of the flex-arm axis can be rotated 360 degrees and extended for further reach. The MS53BA2 flex-arm is manufactured from durable, lightweight recycled aluminum and has the ability to be folded flat, which makes the MS53BA2 stand portable and ready for travel.

## Especificación

| Compatibility        |                                   |
|----------------------|-----------------------------------|
| Compatible with:     | All Dino-Lite models              |
| Housing              |                                   |
| Dimensions           | 35.5cm (W) x 16.5cm (L) x 5cm (H) |
| Weight               | 4 kg                              |
| Features             |                                   |
| ESD safe             | No                                |
| Information          |                                   |
| Package contents     | Stand, table base, HL2 holder     |
| Warranty information | 2 years European warranty         |
| Accessory range:     | Professional stands               |
| Price range          | €400,00 - €600,00                 |

## Product Gallery

