

# RK-04



**Marca:** Dino-Lite  
**Código del producto:** RK-04

## Short Description

Portable rack  
Fast assembly

## Descripción

The Dino-Lite RK-04 stand has a robust, lightweight and compact rack design. It is designed with portability and convenience in mind. All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-04F with fine focus adjustment](#)

## Especificación

|                  |   |
|------------------|---|
| Compatibility    |   |
| Compatible with: | All Dino-Lite models, except Long Working Distance models |
| Housing          |   |
| Dimensions       | 14cm (W) x 11cm (L) x 16.5cm (H)                          |
| Weight           | 605 g   |
| Features         |   |

|                      |                           |
|----------------------|---------------------------|
| ESD safe             | No                        |
| <b>Information</b>   |                           |
| Package contents     | Stand, HD-P1 holder       |
| Warranty information | 2 years European warranty |
| Accessory range:     | Professional stands       |
| Price range          | €100,00 - €200,00         |

### Product Gallery

