

MS53B

Merk: Dino-Lite

Productcode: MS53B

Korte omschrijving

Flexible arm stand

Omschrijving

Dino-Lite heavy duty jointed flex arm stand. The MS53B is a jointed articulating flex-arm type stand with holster for use with the Dino-Lite hand held digital microscopes.

The precision machine crafted components provide smooth adjustment and unlimited horizontal positioning in a well balanced ergonomic design. The MS53B can be easily clamped onto and secured to any workstation or table top. Manufactured from durable, lightweight recycled aluminum and able to be folded flat makes the versatile MS53B stand conveniently portable.

The available MS53BA2 option includes the MS53B flex-arm stand with an additional rigid floor stand and base assembly for use when a workstation table top is not convenient or available.

Specificatie

Compatibiliteit	
Compatibel met:	All Dino-Lite models
Behuizing	
Afmetingen	35.5cm (W) x 16.5cm (L) x 5cm (H)
Gewicht	1.8 kg
Kenmerken	
ESD veilig	No
Informatie	
Inhoud van de verpakking	Stand, HL2 holder
Garantie-informatie	2 years European warranty
accessoires:	Professional stands
Prijsklasse	€350,00 - €500,00

Product galerij

