

# RK-05



**Merk:** Dino-Lite  
**Productcode:** RK-05

## Korte omschrijving

Portable rack  
Tilt and stay feature (30 degrees)  
Fast assembly

## Omschrijving

If you want to observe your object or specimen from a different perspective, the Dino-Lite RK-05 stand provides a multi-angle view with its 30 degrees inclinable pole. RK-05 is a light-weight, ultra-portable stand with a one-step assemble/disassemble design.

All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-05F with fine focus adjustment](#)

## Specificatie

Compatibiliteit	
Compatibel met:	All Dino-Lite models, except Long Working Distance models
Behuizing	

Afmetingen	14cm (W) x 11cm (L) x 16.5cm (H)
Gewicht	611 g
Kenmerken	
ESD veilig	No
Informatie	
Inhoud van de verpakking	Stand, HD-P1 holder
Garantie-informatie	2 years European warranty
accessoires:	Professional stands
Prijsklasse	€150,00 - €250,00

## Product galerij

