

# RK-05



**Merk:** Dino-Lite  
**Productcode:** RK-05

## Korte omschrijving

Portable rack  
Tilt and stay feature (30 degrees)  
Fast assembly

## Omschrijving

If you want to observe your object or specimen from a different perspective, the Dino-Lite RK-05 stand provides a multi-angle view with its 30 degrees inclinable pole. RK-05 is a light-weight, ultra-portable stand with a one-step assemble/disassemble design.

All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-05F with fine focus adjustment](#)

## Specificatie

|                 |   |
|-----------------|---|
| Compatibiliteit |   |
| Compatibel met: | All Dino-Lite models, except Long Working Distance models |
| Behuizing       |   |

|                          |                                  |
|--------------------------|----------------------------------|
| Afmetingen               | 14cm (W) x 11cm (L) x 16.5cm (H) |
| Gewicht                  | 611 g                            |
| Kenmerken                |                                  |
| ESD veilig               | No                               |
| Informatie               |                                  |
| Inhoud van de verpakking | Stand, HD-P1 holder              |
| Garantie-informatie      | 2 years European warranty        |
| accessoires:             | Professional stands              |
| Prijsklasse              | €150,00 - €250,00                |

## Product galerij

