

# RK-06A



**Merk:**  
**Productcode:** RK-06A

## **Korte omschrijving**

Portable rack  
Fast assembly  
Fine-focus adjustment

## **Omschrijving**

The Dino-Lite RK-06A stand is a sturdy and stable mid-range stand solution that can be used with all Dino-Lite digital microscopes. It is constructed of stainless steel and lightweight aluminum and offers precise fine-focus adjustment as well as a quick vertical release function.

Its adjustable horizontal arm has a 15cm range with 360-degree rotation and with the RK-10-EX arm (included) its effective range can be extended to observe larger objects. Various optional accessories are available for the RK-06A, such as the RK-10-FX flexible arm or the RK-10-PX, XY arm. With these features, the RK-06A is the ideal mid-range stand solution for precise orientation and maximum flexibility.

[AVAILABLE EXTENSIONS FOR THE RK-06A STAND](#)

---

**RK-10-EX**



Additional horizontal arm (included with RK-10A)

**RK-10-FX**



Flexible arm extension

**RK-10-VX**



Vertical arm extension

**RK-10-PX2**



X and Y positioning arm

### Specificatie

|                 |                                   |
|-----------------|-----------------------------------|
| Compatibiliteit |                                   |
| Compatibel met: | All Dino-Lite models              |
| Behuizing       |                                   |
| Afmetingen      | 22 cm (L) x 15 cm (W) x 27 cm (H) |
| Gewicht         | 1060 g                            |
| Kenmerken       |                                   |
| ESD veilig      | No                                |

| Informatie               |  |
|--------------------------|--|
| Inhoud van de verpakking | Stand, HD-P1 holder, RK-10-EX horizontal arm |
| Garantie-informatie      | 2 years European warranty                    |
| accessoires:             | Professional stands                          |
| Prijsklasse              | €200,00 - €300,00                            |

### Product galerij

