

SS-M2

Merk:

Productcode: SS-M2



Korte omschrijving

Compatible with RK-10 rack.

Upright or inverted setup

Flexible petri dish holder

Omschrijving

In need of an upright or inverted setup for observing petri dish samples? SS-M2 is a flexible petri dish holder with a two-way configuration setup to facilitate inverted and upright viewing with the use of any Dino-Lite microscope. SS-M2 is only compatible with RK-10 rack.

Specificatie

Behuizing	
Afmetingen	27.5 cm (L) x 15 cm (W) x 4 cm (H)
Gewicht	approx. 240 g = 0.5 lb.

Product galerij

