

# RK-04



**Merk:** Dino-Lite  
**Productcode:** RK-04

## Korte omschrijving

Portable rack  
Fast assembly

## Omschrijving

The Dino-Lite RK-04 stand has a robust, lightweight and compact rack design. It is designed with portability and convenience in mind. All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-04F with fine focus adjustment](#)

## Specificatie

Compatibiliteit	
Compatibel met:	All Dino-Lite models, except Long Working Distance models
Behuizing	
Afmetingen	14cm (W) x 11cm (L) x 16.5cm (H)
Gewicht	605 g
Kenmerken	

ESD veilig	No
<b>Informatie</b>	
Inhoud van de verpakking	Stand, HD-P1 holder
Garantie-informatie	2 years European warranty
accessoires:	Professionele stands
Prijsklasse	€100,00 - €200,00

### Product galerij

