

# RK-04F



**Merk:** Dino-Lite  
**Productcode:** RK-04F

## Korte omschrijving

Portable rack  
Fast assembly  
Fine-focus adjustment

## Omschrijving

The Dino-Lite RK-04F stand has a robust, lightweight and compact rack design. Additionally it has an integrated fine focus adjustment knob to make focusing easier and more precise. The RK-04F is the perfect accessory for users requiring flexibility.

All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-04 without fine focus adjustment](#)

## Specificatie

Compatibiliteit	
Compatibel met:	All Dino-Lite models, except Long Working Distance models
Behuizing	

Afmetingen	14cm (W) x 11cm (L) x 16.5cm (H)
Gewicht	630 g
<b>Kenmerken</b>	
ESD veilig	No
<b>Informatie</b>	
Inhoud van de verpakking	Stand, HD-P1 holder
Garantie-informatie	2 years European warranty
accessoires:	Professional stands
Prijsklasse	€150,00 - €250,00

## Product galerij

