


# CS-41

**Merk:** Dino-Lite  
**Productcode:** CS-41



## Korte omschrijving

Glass calibration slide

## Omschrijving

The CS-41 is a calibration slide that can be used to calibrate your Dino-Lite Edge digital microscope or to be used with the Dino-Eye Edge eyepiece cameras. The total length is 10mm (10.000um) and the smallest gradation on the slide is 0.1mm (100um). The CS-41 calibration slide also contains QR codes in 3 different sizes that can be used for the Auto Calibration procedure within DinoCapture 2.0 software version 1.5.25 or newer. Auto Calibration provides an easy way to calibrate a Dino-Lite Edge serie microscope. The DinoCapture software identifies the special QR codes and performs the calibration automatically.

Row A: magnifications of 120x and above

Row B: magnifications of 40x - 220x

Row C: magnifications of 20x - 50x

[Click here to download the Auto Calibration procedure](#)

## Specificatie

Compatibiliteit	
Compatibel met:	All Dino-Lite models

<b>Behuizing</b>	
Afmetingen	75mm (L) x 25mm (W) x 1,5mm (H)
Gewicht	10 g
<b>Kenmerken</b>	
ESD veilig	No
<b>Informatie</b>	
Inhoud van de verpakking	1 piece
Garantie-informatie	2 years European warranty
accessoires:	Various
Prijsklasse	€0 - €100,00