

AM4023

Merk: Dino-Lite

Productcode: AM4023



Korte omschrijving

1.3 Megapixel resolution
Fits 23mm eyepiece tubes
Advanced software
Including measurements

Omschrijving

The Dino-Eye is a range of digital eyepiece cameras. You can replace the ocular (eyepiece) of regular microscopes with the Dino-Eye which turns the regular microscope into a digital microscope. The AM4023 connects via USB to your PC and fits 23mm eyepiece tubes. The DinoCapture software is bundled with the product.



USB 2.0



Measurement
functionality



1,3 Megapixels
1280 x 1024



Fits 23 mm
oculars

Specificatie

Optiek	
Magnification	Yes, built-in
Type lens	Glass with Anti-Reflection coating

Sensor	
Sensortype	CMOS
Resolutie	1.3 Megapixel (1280x1024)
Maximale framesnelheid	30 fps
Compatibiliteit	
Interface	USB 2.0
Besturingssysteem	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Ondersteunde beeldformaten (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Ondersteunde videoformaten (Windows)	WMV, FLV, SWF
Ondersteunde beeldformaten (MacOS)	JPEG, PNG
Ondersteunde videoformaten (MacOS)	MOV
Beeldvormingsnormen	DirectShow, UVC
Wifi	No
Behuizing	
Materiaal van de behuizing	Metal body
Afmetingen	48 mm (H) x 23 mm (D)
Gewicht	95 g
Kabellengte	1.8 m
Mounting type:	Inside ocular/tube
Inner diameter	23 mm
Kenmerken	
Meting	Yes
Kalibratie	Yes
Microtouch-sensor	No
Informatie	
Inhoud van de verpakking	DinoEye eyepiece camera, carry pouch, software CD, calibration target, user manual
Garantie-informatie	2 years European warranty
Regelgevende goedkeuring	CE, FCC, ROHS
accessoires:	DinoEye Eyepieces
Prijsklasse	€150,00 - €250,00

Product galerij

