

AM4117MTW-FJ



Merk:
Productcode: AM4117MTW-FJ

Korte omschrijving

1.3 Megapixel
High Sensitivity Sensor
Macro Zoom Optics
Adjustable LED Intensity
Robust Metal Housing

Omschrijving

Equipped with 940 nm LEDs, AM4117MTW-FJ suits the examination of surfaces exhibiting distinct optical properties under infrared light. This microscope provides macro-level imaging with zoom-in capabilities, allowing for detailed observation of a wide range of subjects.



1,3 Megapixels
1280 x 960



USB 2.0

Working distance/field of view/depth of field

T Zone

M	WD	FOV (x)	FOV (y)	DOF
50	21	8	6	3
40	42	10	7	6
30	71	13	10	11
20	127	20	15	26
15	183	26	20	47

W Zone

M	WD	FOV (x)	FOV (y)	DOF
30	20	13	10	7
20	33	20	15	14
15	48	26	20	24
10	77	39	29	54
5	166	78	59	266

Specificatie

Verlichting	
Licht / LED-type	940 nm LED
Optiek	
Magnification	10x ~ 50x
Sensor	
Resolutie	1.3 Megapixel (1280 x 960)

Maximale framesnelheid	30 fps at 1.3 MP, MJPEG 30 fps at VGA, MJPEG
Compatibiliteit	
Interface	USB 2.0
Besturingssysteem	Windows 10/11 macOS 10.14 or later
Behuizing	
Vergrotingslot	Yes
Afmetingen	104 mm (H) x 32 mm (D)
Gewicht	Approx. 89 g
Kabellengte	Approx. 1.8 m
Kenmerken	
Meting	Yes
Kalibratie	Yes
Microtouch-sensor	Yes
Informatie	
Inhoud van de verpakking	Microscope, carry pouch, calibration target, user manual, N3C-O- Open cap, N3C-C- Closed cap, N3C-D- diffuser cap, N3C-E- Extension cap, N3C-L- Long cap, N3C-S- Side light cap
Garantie-informatie	2 years European warranty
Regelgevende goedkeuring	CE, FCC, ROHS
Prijsklasse	€800-1000