

# IriScope (MEDL4R)



**Marka:** Dino-Lite Medical  
**Ürün Kodu:** MEDL4R

## Short Description

1.3 Megapixel resolution

Make detailed pictures of the Iris

10-20x Magnification


Medical Device Class 1

## Açıklama

The IriScope (MEDL4R) offers a new, fast and easy way to make pictures of eyes. The IriScope has 2x2 LED's (white/yellow) that illuminate the side of the eye which reduces the reflections in the image as well as the discomfort by bright light. The yellow LED's provide better images of dark irises. The images have a maximum resolution of 1.3 Megapixel.



USB 2.0



Switchable  
white/yellow  
LEDs



Adjustable  
~10-20x



Magnification  
lock



1.3 megapixel  
1280 x 1024

## Özellikler

Lighting

Light/ LED type	White/ Yellow
Number of LEDs	2+2
LED on/off switchable:	Yes
Infrared filter	IR cut-filter >650 nm
Diffuser available	No
Polarizer	No
<b>Optics</b>	
Magnification	10-20x
Lens type	Glass with anti-reflection coating
<b>Sensor</b>	
Sensor type	CMOS
Resolution	1.3 Megapixel (1280x1024)
Maximum frame rate	30 fps
<b>Compatibility</b>	
Interface	USB 2.0
Operating system	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Supported image formats (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Supported video formats (Windows)	WMV, FLV, SWF
Supported image formats (MacOS)	JPEG, PNG
Supported video formats (MacOS)	MOV
Imaging standards	DirectShow, UVC
<b>Housing</b>	
Housing material	Composite/plastic housing
Magnification lock	Yes
Dimensions	10.5cm (L) x 3.2cm (D)
Weight	95 g
Cable length	1.8m
<b>Features</b>	
Special feature	No
Measurement	No
Calibration	No
Microtouch sensor	Yes
<b>Information</b>	
Package contents	Microscope, carry pouch, software CD, user manual
Warranty information	2 years European warranty

Regulatory approval	Medical Device Class 1 – Medical Devices Regulation (EU) 2017/745
Price range	€500,00 - €700,00

**Product Gallery**

