

# RK-04F



**Marka:** Dino-Lite  
**Ürün Kodu:** RK-04F

## Short Description

Portable rack  
Fast assembly  
Fine-focus adjustment

## Açıklama

The Dino-Lite RK-04F stand has a robust, lightweight and compact rack design. Additionally it has an integrated fine focus adjustment knob to make focusing easier and more precise. The RK-04F is the perfect accessory for users requiring flexibility.

All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-04 without fine focus adjustment](#)

## Özellikler

|                  |   |
|------------------|---|
| Compatibility    |   |
| Compatible with: | All Dino-Lite models, except Long Working Distance models |
| Housing          |   |

|                      |                                  |
|----------------------|----------------------------------|
| Dimensions           | 14cm (W) x 11cm (L) x 16.5cm (H) |
| Weight               | 630 g                            |
| Features             |                                  |
| ESD safe             | No                               |
| Information          |                                  |
| Package contents     | Stand, HD-P1 holder              |
| Warranty information | 2 years European warranty        |
| Accessory range:     | Professional stands              |
| Price range          | €150,00 - €250,00                |

## Product Gallery

